

# Endodontic Focus

Featuring products from SybronEndo, Coltene, Quala and more!

## K3<sup>TM</sup> XF NiTi Rotary Files

Buy 50 packs, get 25 packs **FREE!**

WE SUBMIT THE COUPON FOR YOUR **FREE GOODS!**

As low as **55.95** (25+) per 6-pack!

See inside page for details.

**SybronEndo**  
Sybron Dental Specialties

\*Offer valid through March 31, 2012.



## K3<sup>®</sup> NiTi Rotary Files



As low as **55.95** (25+) per 6-pack!


See inside page for details.

**SybronEndo**  
Sybron Dental Specialties

## HYFLEX<sup>®</sup> CM NiTi Rotary Files

WE SUBMIT THE COUPON FOR YOUR **FREE GOODS!**

Buy 30 refills, get 8 same **FREE!**



Only **55.95** per 6-pack!  
See back cover for details.

**COLTENE<sup>®</sup>**

## AH Plus<sup>®</sup> Root Canal Sealer

Buy 2 kits, get an AH Plus Jet starter kit **FREE!**

WE SUBMIT THE COUPON FOR YOUR **FREE GOODS!**



(160-7020) Kit **107.95**  
Contains: 4ml paste A, 4ml paste B, mixing pad and directions for use.

**DENTSPLY**  
MAILLEFER

## Quala Tapered Gutta Percha Points



Only **7.50** (6+) per 6-pack!

See back cover for details.

**Quala**  
DENTAL PRODUCTS

## 17% EDTA Solution (Pulpdent)



Only **19.50** per 120ml bottle!  
(570-1860)

## Glyde-File Prep (Dentsply Maillefer)



(570-1800) Syringe Kit **30.50**  
Contains: (5) 3ml syringes and 50 disposable tips.

# MidwestDental<sup>®</sup>

EQUIPMENT • SERVICES • SUPPLIES

Phone (800) 766-2025 | [www.midwestdental.com](http://www.midwestdental.com)



# K3<sup>®</sup> Nickel Titanium Rotary Files

Each Assorted Pack . . . . . **60.95**  
 25+ . . . . . **55.95**

## K3 Procedure Pack

Contains: One each of .10 taper 17mm #25, .08 taper 17mm #25, 25, 30, 35 and 40.

(332-4450) .04 21mm      (332-4455) .04 25mm  
 (332-4460) .06 21mm      (332-4465) .06 25mm

## K3 VTVT Pack

Contains: One each of .10 taper #25, .08 taper #25, .06 taper #35, .06 taper #25, .04 taper #30 and .04 taper #20.

(332-4485) 21mm      (332-4490) 25mm

## K3 G-Pack

Contains: One each of .12 taper 17mm #25, .10 taper #25, .08 taper #25, .06 taper #25, .04 taper #25 and .02 taper #25.

(332-4510) 21mm      (332-4515) 25mm

Taper	Refill			
	Tip	17mm	21mm	25mm
.08	25	332-4505	332-4495	332-4500
.10	25	332-4000	332-4440	332-4445
.12	25	332-4480	332-4470	332-4475
<b>Pricing</b>	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95</b>			

Tip	.04 Taper Refill		
	21mm	25mm	30mm
15	332-4135	332-4140	332-4130
20	332-4150	332-4155	332-4145
25	332-4165	332-4170	332-4160
30	332-4180	332-4185	332-4175
35	332-4195	332-4200	332-4190
40	332-4210	332-4215	332-4205
45	332-4225	332-4230	332-4220
50	332-4240	332-4245	332-4235
55	332-4255	332-4260	332-4250
60	332-4270	332-4275	332-4265
15-40	332-4120	332-4125	
<b>Pricing</b>	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95</b>		



K3 is a triple-fluted, asymmetric endo file system designed to cut quickly, efficiently and safely, with unparalleled debris removal. And its unique, cross-sectional design makes K3 files significantly more resistant to cyclic fatigue than some other rotary files. K3 files are able to negotiate even the most challenging S-shape curves. They do so with little or no canal deformation. And the safe-ended tip assists in following the canal shape while minimizing apical transportation.

K3 files are available in 6 different tapers with a multitude of lengths and tip sizes. This allows you to customize your technique for any case from single, straight roots to calcified, curved molars. K3 files are also available in a number of pre-packaged assorted packs, which is a perfect way to get started.

Tip	.02 Taper Refill		
	21mm	25mm	30mm
15	332-4020	332-4025	332-4015
20	332-4035	332-4040	332-4030
25	332-4050	332-4055	332-4045
30	332-4065	332-4070	332-4060
35	332-4080	332-4085	332-4075
40	332-4095	332-4100	332-4090
45	332-4110	332-4115	332-4105
15-40	332-4005	332-4010	
<b>Pricing</b>	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95</b>		

Tip	.06 Taper Refill		
	21mm	25mm	30mm
15	332-4295	332-4300	332-4290
20	332-4310	332-4315	332-4305
25	332-4325	332-4330	332-4320
30	332-4340	332-4345	332-4335
35	332-6351	332-4360	332-4350
40	332-4370	332-4375	332-4365
45	332-4385	332-4390	332-4380
50	332-4400	332-6505	332-4395
55	332-4415	332-4420	332-4410
60	332-4430	332-4435	332-4425
15-40	332-4280	332-4285	
<b>Pricing</b>	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95</b>		



## K3<sup>TM</sup> XF Nickel Titanium Rotary Files

K3XF provides the safety and self-centering features of the original K3, plus an extraordinary new level of flexibility and resistance to cyclic fatigue provided by SybronEndo's proprietary R-Phase<sup>TM</sup> Technology. K3XF is the ultimate in safety, flexibility and control for rotary file users looking for a stronger and more flexible file.

K3XF features an Axxess stainless steel handle for more working room, third radial land for centered stability, and a safe-ended tip assists in following the canal shape while minimizing apical transportation. The file's variable pitch also provides the ultimate in control over other NiTi files.

### K3XF is engineered for efficiency:

- Positive rake angle provides the active cutting action of the K3 file
- The third radial land stabilizes and keeps the instrument centered in the canal and minimizes "over-engagement"
- Reduced radial land with blade support adds peripheral strength to resist torsional and rotary stresses
- Radial land relief reduces friction on the canal wall

### X3XF Files feature the following benefits:

- Heat treated with SybronEndo's proprietary R-Phase heat treatment process for half the stiffness of other files and two times more resistant to fracture
- Third radial land for better centering, more stability and less chance of transportation
- Variable pitch for more control and assists the file in channeling debris coronally
- Axxess stainless steel handle for better posterior access



## SPECIAL OFFERS!

Buy 50 packs, get 25 packs FREE!<sup>\*</sup>  
**OR**  
 Buy 15 packs, get 4 packs FREE,  
**PLUS 2 packs of K3**  
**Gutta Percha Points FREE!**

Call us at (800) 766-2025  
 or contact your dental products  
 specialist for additional offers!



Taper	Refill			
	Tip	17mm	21mm	25mm
.08	25	332-5295	332-5285	332-5290
.10	25	332-5010	332-5000	332-5005
Pricing	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95*</b> Net <b>37.30*</b>			

Tip	.04 Taper Refill		
	21mm	25mm	30mm
20	332-5020	332-5025	332-5015
25	332-4251	332-5040	332-5030
30	332-4301	332-5055	332-5045
35	332-5065	332-4355	332-5060
40	332-5080	332-4405	332-5075
45	332-5095	332-5100	332-5090
50	332-5110	332-5115	332-5105
55	332-5125	332-5130	332-5120
60	332-5140	332-5145	332-5135
Pricing	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95*</b> Net <b>37.30*</b>		

Tip	.06 Taper Refill		
	21mm	25mm	30mm
20	332-5155	322-6205	332-5150
25	332-5170	332-5175	332-5165
30	322-5185	322-6305	332-5180
35	332-5200	322-6355	332-5195
40	332-5215	332-5220	332-5210
45	332-5230	332-5235	332-5225
50	332-5245	332-5250	332-5240
55	332-5260	332-5265	332-5255
60	332-5275	332-5280	332-5270
Pricing	Pkg/6 ... <b>60.95</b> , 25+ <b>55.95*</b> Net <b>37.30*</b>		

\*50+25 offer valid through March 31, 2012.

# MidwestDental®

EQUIPMENT • SERVICES • SUPPLIES

**COLTENE®**



**SPECIAL OFFER!**  
Buy 30 refills, get 8 same FREE!

## **HYFLEX® CM** Controlled Memory Nickel Titanium Files



(332-3500) Intro Kit . . . . . **347.95**

Contains: Six rotary file assorted packs (2-21mm, 4-25mm), Spectra Point® GT Gutta Percha Points, one endo procedure block, three endo file practice blocks, CD and a step-by-step card.

(332-3665) Assortment . . . . . **55.95**

Contains: One each of 04/20 21mm, 04/25 21mm, 04/30 21mm, 04/40 21mm, 06/20 21mm and 08/25 19mm.

(332-3670) Assortment . . . . . **55.95**

Contains: One each of 04/20 25mm, 04/25 25mm, 04/30 25mm, 04/40 25mm, 06/20 25mm and 08/25 19mm.

(332-3505) 08/25 19mm, pkg/6. . . . . **55.95**

(332-3655) Practice Blocks, pkg/3. . . . . **22.95**

(332-3660) Procedure Block, each. . . . . **41.50**

Tip	.04 Taper Refill	
	21mm	25mm
15	332-3510	332-3585
20	332-3515	332-3590
25	332-3520	332-3595
30	332-3525	332-3600
35	332-3530	332-3605
40	332-3535	332-3610
45	332-3540	332-3615
50	332-3545	332-3620
60	332-3550	332-3625
<b>Pricing</b>	Pkg/6 ... <b>55.95</b> , Net <b>44.17</b>	

Tip	.06 Taper Refill	
	21mm	25mm
20	332-3560	332-3630
25	332-3565	332-3635
30	332-3570	332-3640
35	332-3575	332-3645
40	332-3580	332-3650
<b>Pricing</b>	Pkg/6 ... <b>55.95</b> , Net <b>44.17</b>	



## **Quala Tapered Gutta Percha Points**

For use with matching .04 and .06 tapered files in conjunction with single-cone obturating techniques. Color-coded ISO sizes for fast, easy identification. Hand rolled.



Standardized ISO Sizes .04 Taper  
Individual sizes, box/100 . . . . . **7.95**  
6+ . . . . . **7.50**

(338-3600) 15	(338-3620) 35	(338-3640) 55
(338-3605) 20	(338-3625) 40	(338-3645) 60
(338-3610) 25	(338-3630) 45	(338-3645) 15-40
(338-3615) 30	(338-3635) 50	

Standardized ISO Sizes .06 Taper  
Individual sizes, box/60 . . . . . **7.95**  
6+ . . . . . **7.50**

(338-3655) 15	(338-3670) 30	(338-3685) 45
(338-3660) 20	(338-3675) 35	(338-3690) 50
(338-3665) 25	(338-3680) 40	(338-3695) 15-40

**Place your order today!**

**Phone (800) 766-2025 | www.midwestdental.com**

Prices valid January 1 through June 30, 2012. Not responsible for printing and typographical errors. All promotions are subject to current manufacturer specials. Code: MWDENDO12